



# Recombinant Human MUC-1 (a.a 24-380) [hFc] (DAG-WT2244)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	hFc
Applications	Immunoassays
Molecular Weight	60.6 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Introduction	Mucin-1 (MUC-1) is a heterodimer transmembrane protein of the mucin family encoded in humans by the MUC1 gene. It is cleaved into two chains: mucin-1 subunit alpha (MUC1-NT; MUC1-alpha) and mucin-1 subunit beta (MUC-CT; MUC1-beta). These subunits differ in size due to proteolytic cleavage of the translated precursor protein in the endoplasmic reticulum.
Keywords	MUC1; Mucin-1

# GENE INFORMATION

UniProt ID	P15941
------------	--------

---