



Recombinant Human IGFBP-5 [His] (DAG-WT2247)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	29.6 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM PB, 150mM NaCl. pH7.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.
Keywords	Insulin-like growth factor-binding protein 5; IGFBP5

GENE INFORMATION

UniProt ID

P24593
