



Recombinant Human PRTN3 [His] (DAG-WT2251)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	25.64 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	10mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Proteinase 3, also known as PRTN3, is an enzyme that in humans is encoded by the PRTN3 gene.
Keywords	Myeloblastin; PRTN3

GENE INFORMATION

