



Recombinant Human GM-CSF (DAG-WT2255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	14.5 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	10mM Tris-HCl, 4% Mannitol, 1% Sucrose, pH 8.5
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes
Keywords	Granulocyte-macrophage colony-stimulating factor; CSF2

GENE INFORMATION

UniProt ID	P04141
------------	--------