



Recombinant Human Amyloid β A4 Protein [His] (DAG-WT2263)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	33.6 kDa
Format	Liquid
Size	50 μ g, 500 μ g, 1 mg
Buffer	20mM PB, 150mM NaCl, pH7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Functions as a cell surface receptor and performs physiological functions on the surface of neurons relevant to neurite growth, neuronal adhesion and axonogenesis. Interaction between APP molecules on neighboring cells promotes synaptogenesis.
Keywords	Amyloid-beta precursor protein; APP

GENE INFORMATION

