



Recombinant Human Vitronectin [His] (DAG-WT2265)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	53.35 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM Tris-HCl, 150mM NaCl, pH 8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Vitronectin is a multifunctional glycoprotein present in blood and in the extracellular matrix. It binds glycosaminoglycans, collagen, plasminogen and the urokinase-receptor, and also stabilizes the inhibitory conformation of plasminogen activation inhibitor-1.
Keywords	Vitronectin; VN

GENE INFORMATION

