



Geneticin [BSA] (DAG-WT2284)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Cas.No	108321-42-2
Applications	ELISA, LF
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	G418 is an aminoglycoside antibiotic similar in structure to gentamicin B1. It is produced by <i>Micromonospora rhodorangea</i> . G418 blocks polypeptide synthesis by inhibiting the elongation step in both prokaryotic and eukaryotic cells.
Keywords	Geneticin; Aminoglycoside antibiotic; <i>Micromonospora rhodorangea</i> ; G418