



Puromycin [BSA] (DAG-WT2285)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LF
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Puromycin is an antibiotic protein synthesis inhibitor which causes premature chain termination during translation.
Keywords	Puromycin; Aminonucleoside antibiotic; Streptomyces alboniger