



Recombinant BHV-1 Envelope Glycoprotein E Protein [His] (DAG-WT236)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant BHV-1 Envelope glycoprotein E is produced by Mammalian expression system and expressed with a 6His tag at the C-terminus.
Purity	> 95% , as determined by SDS-PAGE
Conjugate	His
Applications	Immunogen, calibrator / Standard
Molecular Weight	43kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	PBS, PH 7.4
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Introduction	Bovine alphaherpesvirus 1 (BoHV-1) is a virus of the family Herpesviridae and the subfamily Alphaherpesvirinae, known to cause several diseases worldwide in cattle, including rhinotracheitis, vaginitis, balanoposthitis, abortion, conjunctivitis, and enteritis.
Keywords	Bovine alphaherpesvirus 1; BHV-1; BoHV-1; Herpesviridae