



Glycylproline dipeptidyl aminopeptidase (GPDA) control (DAG-WT2467)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥250U/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Glycylproline dipeptidyl aminopeptidase; GPDA
----------	---