



## Isozyme LDH-1 control (DAG-WT2475)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                       |
|---------------|-----------------------|
| Conjugate     | Unconjugated          |
| Applications  | Immunoassays          |
| Format        | Liquid or Lyophilized |
| Concentration | ≥300U/L               |
| Size          | 0.5 mL, 1.0 mL        |
| Buffer        | Synthetic serum       |
| Storage       | Store at -20°C        |

### BACKGROUND

|          |   |
|----------|---|
| Keywords | Isozyme LDH-1; Lactate dehydrogenase; LDH |
|----------|---|