



α -Hydroxybutyrate dehydrogenase (HBDH) control (DAG-WT2476)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------|--------------------------------|
| Conjugate | Unconjugated |
| Applications | Immunoassays |
| Format | Liquid or Lyophilized |
| Concentration | $\geq 500\text{U/L}$ |
| Size | 0.5 mL, 1.0 mL |
| Buffer | Synthetic serum |
| Storage | Store at -20°C |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Alpha-hydroxybutyrate dehydrogenase is a marker of cell death reflecting red blood cell, myocardial, liver and renal damage. It occurs in higher levels in heart muscle and can be used to aid the diagnosis of myocardial infarction |
| Keywords | α -Hydroxybutyrate dehydrogenase; HBDH |