



## Creatine kinase (CK) control (DAG-WT2478)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥1000U/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

### BACKGROUND

Introduction	Creatine kinase-MB (CK-MB) is a form of an enzyme found primarily in heart muscle cells. This test measures CK-MB in the blood. CK-MB is one of three forms (isoenzymes) of the enzyme creatine kinase (CK). These isoenzymes include: CK-MM (found in skeletal muscles and the heart)
Keywords	Creatine kinase MB; CKMB