



Lp-PLA2 control (DAG-WT2497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥350ng/mL
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Lipoprotein-associated phospholipase A2L Lp-PLA2; Platelet-activating factor acetylhydrolaseL PAF-AH
----------	---