



# Glucose-6-phosphate dehydrogenase (G6PD) control (DAG-WT2498)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥2000U/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

## BACKGROUND

Keywords	Glucose-6-phosphate dehydrogenase; (G6PD)
----------	---