



# Thymidine kinase 1 (TK1) control (DAG-WT2501)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                       |
|---------------|-----------------------|
| Conjugate     | Unconjugated          |
| Applications  | Immunoassays          |
| Format        | Liquid or Lyophilized |
| Concentration | 15pmol/L              |
| Size          | 0.5 mL, 1.0 mL        |
| Buffer        | Synthetic serum       |
| Storage       | Store at -20°C        |

## BACKGROUND

|          |                         |
|----------|-------------------------|
| Keywords | Thymidine kinase 1; TK1 |
|----------|-------------------------|