



## Indirect bilirubin control (DAG-WT2508)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                       |
|---------------|-----------------------|
| Conjugate     | Unconjugated          |
| Applications  | Immunoassays          |
| Format        | Liquid or Lyophilized |
| Concentration | ≥400μmol/L            |
| Size          | 0.5 mL, 1.0 mL        |
| Buffer        | Synthetic serum       |
| Storage       | Store at -20°C        |

### BACKGROUND

|          |               |
|----------|---------------|
| Keywords | Bilirubin; BR |
|----------|---------------|