



# Triglyceride (TAG) control (DAG-WT2511)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥20mmol/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

## BACKGROUND

Introduction	<p>A triglyceride is an ester derived from glycerol and three fatty acids (from tri- and glyceride). Triglycerides are the main constituents of body fat in humans and other vertebrates, as well as vegetable fat. They are also present in the blood to enable the bidirectional transference of adipose fat and blood glucose from the liver, and are a major component of human skin oils.</p>
Keywords	Triglyceride; TAG