



Apolipoprotein AII (ApoA-II) control (DAG-WT2515)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	1.0g/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism
Keywords	Apolipoprotein A-II; APOA2