



# Apolipoprotein C-III (ApoC-III) control (DAG-WT2517)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	400mg/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

## BACKGROUND

Keywords	Apolipoprotein C-III; ApoC-III
----------	--------------------------------