



Apolipoprotein E (ApoE) control (DAG-WT2521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥150mg/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	Apolipoprotein E (APOE) is a class of apolipoprotein found in the chylomicron and IDLs that binds to a specific receptor on liver cells and peripheral cells. It is essential for the normal catabolism of triglyceride-rich lipoprotein constituents.
Keywords	Apolipoprotein E; APOE