



# Remnant-like particle cholesterol (RLP-C) control (DAG-WT2524)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥30mg/dL (2.5mmol/L)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

## BACKGROUND

Keywords	Remnant-like particle cholesterol; RLP
----------	--