



## Uric acid (UA) control (DAG-WT2527)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                |
|---------------|--------------------------------|
| Conjugate     | Unconjugated                   |
| Applications  | Immunoassays                   |
| Format        | Liquid or Lyophilized          |
| Concentration | $\geq 1200 \mu\text{mol/L}$    |
| Size          | 0.5 mL, 1.0 mL                 |
| Buffer        | Synthetic serum                |
| Storage       | Store at $-20^{\circ}\text{C}$ |

### BACKGROUND

|          |               |
|----------|---------------|
| Keywords | Uric acid; UA |
|----------|---------------|