



Glycated albumin (GA) control (DAG-WT2540)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥250μmol/L (20g/L)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	Human albumin is a protein found in the blood plasma of humans and other mammals. It is produced by liver cells and helps in transportation of drugs, hormones, and essential nutrients throughout the body. Albumin also helps maintain blood volume and pressure. It can be used in medical treatments, such as in the treatment of shock or as a blood substitute. It is also used in cosmetic and beauty treatments.
Keywords	Albumin; ALB