



# Alpha Fetoprotein (AFP) control (DAG-WT2580)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥400IU/mL (500ng/mL)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

## BACKGROUND

Keywords	Alpha Fetoprotein; AFP
----------	------------------------