



## Calcitonin (CT) control (DAG-WT2590)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	5ng/mL
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

### BACKGROUND

Introduction	Calcitonin gene-related peptide (CGRP) is a member of the calcitonin family of peptides consisting of calcitonin, amylin, adrenomedullin, adrenomedullin 2 (intermedin) and calcitonin -receptor-stimulating peptide. Calcitonin is mainly produced by thyroid C cells whilst CGRP is secreted and stored in the nervous system.
Keywords	Calcitonin; Blood calcium; CALCA