



Neutrophil gelatinase-associated lipocalin (NGAL) control (DAG-WT2598)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥2000µg/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	LCN2 (lipocalin 2) is a protein-coding gene. Diseases associated with LCN2 include appendicitis, and proctitis. GO annotations related to this gene include iron ion binding and protein homodimerization activity. An important paralog of this gene is LCN12.
Keywords	LCN2; Lipocalin 2; NGAL