



## Fructosamine (FMN) control (DAG-WT2603)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥4mmol/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

### BACKGROUND

Introduction	Fructosamine (1-amino-1-deoxy fructose), is a stable ketoamine, formed by the reaction between glucose and the amino group of protein (predominantly albumin, but also including globulins and lipoprotein)
Keywords	Fructosamine; FMN