



Alpha-2-macroglobulin (α2MG) control (DAG-WT2607)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	6.5g/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	α2-Macroglobulin; α2M; alpha-2-macroglobulin
----------	--