



Haptoglobin (HPT) control (DAG-WT2609)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	3.5g/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	<p>Haptoglobin is a protein made by your liver. It attaches to a certain type of hemoglobin.</p> <p>Hemoglobin is a protein in your red blood cells that carries oxygen from your lungs to the rest of your body. Most hemoglobin is located inside red blood cells, but small amounts circulate in the bloodstream.</p>
Keywords	Haptoglobin; HPT