



Insulin-like growth factor-1 (IGF-1) control (DAG-WT2623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	100ng/mL
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Insulin-like growth factor-1; IGF-1
----------	-------------------------------------