



1,5-Anhydroglucitol (1,5-AG) control (DAG-WT2629)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	$\geq 300 \mu\text{mol/L}$
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	1,5-Anhydroglucitol; 1,5-AG
----------	-----------------------------