



Recombinant Malus Domestica Mal d 4 [His] (DAG-WT264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant Mal d 4 is produced by E.coli and the target gene encoding Ser2-Leu131 is expressed with a His tag at the N-terminus.
Species	Malus domestica
Purity	> 95% , as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	14 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	50 mM Na ₂ HPO ₄ (pH 7.4), 0.5 M NaCl, 30 mM imidazole, 6 M urea
Preservative	None
Storage	Short Term: 4°C. Long Term: -20°C to -80°C. Avoid repeated freezing and thawing.

BACKGROUND

Keywords	Malus domestica; Mal d 4; Profilin; Mal d 4.0202
----------	--

GENE INFORMATION

UniProt ID	Q84RR6
------------	--------
