



CA-125 control (DAG-WT2642)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥400IU/mL
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	CA-125 has found application as a tumor marker or biomarker that may be elevated in the blood of some patients with specific types of cancers, or other benign conditions. MUC16 (CA125) is a serum marker that is used routinely in gynecologic practice to monitor patients with ovarian cancer.
Keywords	CA125; MUC16; Tumor marker