



Hepatocyte Growth Factor (HGF) control (DAG-WT2646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|-----------------------|
| Conjugate | Unconjugated |
| Applications | Immunoassays |
| Format | Liquid or Lyophilized |
| Concentration | 400ng/L |
| Size | 0.5 mL, 1.0 mL |
| Buffer | Synthetic serum |
| Storage | Store at -20°C |

BACKGROUND

| | |
|----------|-------------------------------|
| Keywords | Hepatocyte Growth Factor; HGF |
|----------|-------------------------------|