



Progesterone control (DAG-WT2658)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥70ng/mL (223nmol/L)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Progesterone; Pregn-4-ene-3,20-dione; P4; Pregnenedione
----------	---