



17-OH-Progesterone control (DAG-WT2659)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥65nmol/L (21.5ng/mL)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	17-hydroxyprogesterone; 17-OHP; 17-OH-Progesterone
----------	--