



Adrenocorticotrophic hormone (ACTH) control (DAG-WT2665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥23pmol/L (102.6pg/mL)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Adrenocorticotrophic hormone; ACTH
----------	------------------------------------