



## Aldosterone (ALD) control (DAG-WT2668)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	$\geq 2\mu\text{mol/mL}$ (720pg/mL)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

### BACKGROUND

Keywords	Aldosterone; ALD
----------	------------------