



# Human collagen type III control (DAG-WT2669)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                       |
|---------------|-----------------------|
| Conjugate     | Unconjugated          |
| Applications  | Immunoassays          |
| Format        | Liquid or Lyophilized |
| Concentration | ≥360ng/mL             |
| Size          | 0.5 mL, 1.0 mL        |
| Buffer        | Synthetic serum       |
| Storage       | Store at -20°C        |

## BACKGROUND

|          |                      |
|----------|----------------------|
| Keywords | COL1A3; Collagen III |
|----------|----------------------|