



# Anti-Müllerian hormone (AMH) control (DAG-WT2674)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥10ng/mL
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

## BACKGROUND

Keywords	Anti-Müllerian hormone; AMH
----------	-----------------------------