



25-hydroxyvitamin D control (DAG-WT2678)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|-----------------------|
| Conjugate | Unconjugated |
| Applications | Immunoassays |
| Format | Liquid or Lyophilized |
| Concentration | ≥240nmol/L |
| Size | 0.5 mL, 1.0 mL |
| Buffer | Synthetic serum |
| Storage | Store at -20°C |

BACKGROUND

| | |
|----------|------------------------------|
| Keywords | 25-hydroxyvitamin D; 25(OH)D |
|----------|------------------------------|