



25-hydroxyvitamin D2 control (DAG-WT2679)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥140nmol/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	25-hydroxyvitamin D2; Ergocalciferol
----------	--------------------------------------