



25-hydroxyvitamin D3 control (DAG-WT2680)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥150nmol/L
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	25-hydroxyvitamin D3; Calcifediol; 25-hydroxycholecalciferol
----------	--------------------------------------------------------------