



## Vitamin B12 control (DAG-WT2686)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥1600pmol/L (2167pg/mL)
Size	0.5 mL
Buffer	Synthetic serum
Storage	Store at -20°C

### BACKGROUND

Introduction	Vitamin B12, also known as cobalamin, is a water-soluble vitamin involved in metabolism. It is one of eight B vitamins. It is required by animals, which use it as a cofactor in DNA synthesis, and in both fatty acid and amino acid metabolism.
Keywords	Vitamin B12; VB12