



Vitamin C control (DAG-WT2687)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	$\geq 150 \mu\text{mol/L}$
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Vitamin C; VC; Ascorbic acid
----------	------------------------------