



Hydrocortisone control (DAG-WT2700)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	≥1800nmol/L (64µg/dL)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Hydrocortisone; A-hydrocort; Cortisol
----------	---------------------------------------