



## Acetaminophen control (DAG-WT2704)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	$\geq 1000 \mu\text{mol/L}$ (150mg/L)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

### BACKGROUND

Keywords	Acetaminophen; para-hydroxyacetanilide
----------	--