



Phenobarbital control (DAG-WT2705)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	$\geq 240 \mu\text{mol/L}$ (55.6 $\mu\text{g/mL}$)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Phenobarbital; Phenobarb
----------	--------------------------