



Salicylic acid control (DAG-WT2706)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|-----------------------|
| Conjugate | Unconjugated |
| Applications | Immunoassays |
| Format | Liquid or Lyophilized |
| Concentration | ≥5mmol/L (69mg/dL) |
| Size | 0.5 mL, 1.0 mL |
| Buffer | Synthetic serum |
| Storage | Store at -20°C |

BACKGROUND

| | |
|----------|---|
| Keywords | Salicylic acid; C ₆ H ₄ (OH)COOH; Aspirin |
|----------|---|